


MOBILE STORAGE SOLUTION FOR INTRAVENOUS POLES



IV POLE STORAGE TROLLEY: IT916-BLS

Intravenous (IV) poles are constantly being removed from beds and require safe storage. Introducing the Polymedic IV Pole Storage Trolley, which allows you to easily store up to 75 IV poles in the one trolley. Developed in conjunction with a major Brisbane hospital, this trolley promotes space saving vertical storage, that allows IV poles to be aired out, vertically stored and conveniently stored without any OH&S issues, offering a safer solution for both staff and patients.

PRODUCT FEATURES:

- Durable, light weight frame for easy use;
- A storage solution that keeps the IV poles off the floor, and mobile;
- Four easy to manoeuvre 75mm castor wheels with braking system;
- Safe storage reduces damage to IV poles, increasing product longevity;
- Numbers of IV poles can be assessed quickly;
- Aeration promotes better product recovery;
- Standard size: 1060(W) x 650(D) x 950(H)mm, housing up to 75 IV poles;
- Trolley can be modified as required;
- OH&S risks reduced;
- Manufactured in Australia.

The Polymedic design team can create products to suit your specific requirements.
For more information, please call us on (07) 5539 3633 or email sales@polymedic.com.au.

Freecall: 1800 426 856 [www: polymedic.com.au](http://www.polymedic.com.au)