


MOBILE STORAGE SOLUTION FOR PATIENT LIFT BARS



PATIENT LIFT BAR STORAGE TROLLEY: IT920-BLS

Patient lift bars are constantly being removed from beds and require safe storage. Introducing the Polymedic Patient Lift Bar Storage Trolley, which allows you to easily store up to 20 bars in the one trolley. Developed in conjunction with a major Brisbane hospital, this trolley promotes space saving vertical storage, that allows patient lift bars to be aired out, vertically stored and conveniently stored without any OH&S issues, offering a safer solution for both staff and patients.

PRODUCT FEATURES:

- Durable, light weight frame for easy use;
- A storage solution that keeps the patient lift bars off the floor, and mobile;
- Four easy to manoeuvre 125mm castor wheels with braking system;
- Safe storage reduces damage to patient lift bars, increasing product longevity;
- Numbers of patient lift bars can be assessed quickly;
- Aeration promotes better product recovery;
- Standard size: 850(W) x 1000(D) x 1305(H)mm, housing up to 20 patient lift bars;
- Trolley can be modified as required;
- OH&S risks reduced;
- Manufactured in Australia.

The Polymedic design team can create products to suit your specific requirements. For more information, please call us on (07) 5539 3633 or email sales@polymedic.com.au.

Freecall: 1800 426 856 [www: polymedic.com.au](http://www.polymedic.com.au)