


MOBILE STORAGE SOLUTION FOR WALKING CRUTCHES



CRUTCHES STORAGE TROLLEY: IT925-BLS

Walking crutches are often left against walls or on the floor, creating unsafe obstacles. The Polymedic Crutches Storage Trolley provides a safe and easy storage solution for walking crutches. Developed in conjunction with a major Brisbane hospital, this trolley promotes space saving vertical storage, that allows walking crutches to be aired out and conveniently stored without any OH&S issues, offering a safer solution for both staff and patients.

PRODUCT FEATURES:

- Durable, light weight frame for easy use;
- A storage solution that keeps the crutches off the floor, and mobile;
- Four easy to manoeuvre 75mm castor wheels with braking system;
- Safe storage reduces damage to crutches increasing product longevity;
- Crutches can be stored and assessed quickly;
- Aeration promotes better product recovery;
- Standard size: 1335(W) x 590(D) x 2100(H)mm, housing up to 30 Bariatric crutches, 40 standard crutches, and 10 arm crutches;
- Trolley can be modified as required;
- OH&S risks reduced;
- Manufactured in Australia.

The Polymedic design team can create products to suit your specific requirements. For more information, please call us on (07) 5539 3633 or email sales@polymedic.com.au.

Freecall: 1800 426 856 [www: polymedic.com.au](http://www.polymedic.com.au)